The new standard in hip sonography Sonoffst

Anatomically formed

fixation areas:

(2) shoulder and arm

(1) hip

(3) head

orthopunkt ag Schöngrünstrasse 35 CH-4500 Solothurn P +41 (0)32 623 00 53 F +41 (0)32 623 00 58 www.orthopunkt.ch

Spring supported guidance of the ultrasonic head

Universal or dedicated holders for ultrasonic head

Controlled guidance secures tilt-free use of the ultrasonic head for repeatable outcome

Disinfectable and urine resisting surface

Sturdy and yet mobile sheet construction

Small foot print

Fast and safe single-handed fixation

Great travel range gives freedom in

all room situations

Modern design and

functional ergonomics

100 % Swiss made

Clean cable runs

Helpful holder for ultrasonic gel supply

Equally usable for left and right handers due to three different installation options

